

# FINNZYMES

# Phusion™

## High-Fidelity PCR Kit

### Product codes:

F-553S, 50 U

F-553L, 200 U

Stable for one year from the packaging date. Store at -20°C.

## 1. Introduction

Finnzymes' Phusion™ High-Fidelity DNA Polymerase offers extreme performance for all PCR applications. Incorporating an exciting new technology, Phusion DNA Polymerase brings together a novel *Pyrococcus*-like enzyme with a processivity-enhancing domain. Phusion DNA Polymerase generates long templates with an accuracy and speed previously unattainable with a single enzyme, even on the most difficult templates. The extreme fidelity makes Phusion DNA Polymerase a superior choice for cloning. Using a *lacI*-based method modified from previous studies,<sup>1</sup> the error rate of Phusion DNA Polymerase in Phusion HF Buffer is determined to be  $4.4 \times 10^{-7}$ , which is approximately 50-fold lower than that of *Thermus aquaticus* DNA polymerase, and 6-fold lower than that of *Pyrococcus furiosus* DNA polymerase.

The Phusion High-Fidelity PCR Kit includes lambda DNA control template and primers for 1.3 kb and 10 kb amplicons. The template amount is sufficient for performing 20 control reactions in a 50 µl reaction volume or 50 control reactions in a 20 µl volume.

**Phusion™ DNA Polymerase is unlike other enzymes. Please read the Quick Guide to modify your protocol for optimal results!**

### Quick Guide:

- Use Phusion DNA Polymerase at 0.5-1.0 U per 50 µl reaction volume. Do not exceed 2 U/50 µl. (see 4.1)
- Use 15-30 s/kb for extension. Do not exceed 1 min/kb. (see 5.4)
- Use 98°C for denaturation. (see 5.1 & 5.2)
- Anneal at  $T_m + 3^\circ\text{C}$  (> 20nt) or use 2-step protocol. Calculate  $T_m$ 's using nearest-neighbor method. (see 5.3)
- Use 200 µM of each dNTP. Do not use dUTP. (see 4.3)
- Note: Phusion DNA Polymerase produces blunt end DNA products.

## 2. Kit Components

Component	Concentration	F-553S	F-553L
Phusion™ DNA Polymerase	2 U/µl	50 U	200 U
5x Phusion™ HF buffer		1 x 1.5 ml	3 x 1.5 ml
5x Phusion™ GC buffer		1 x 1.5 ml	1 x 1.5 ml
dNTP mix	10 mM each	100 µl	400 µl
MgCl <sub>2</sub> solution	50 mM	1.5 ml	1.5 ml
Control lambda template	0.5 ng/µl	40 µl	40 µl
1.3 kb primers	4 µM each	50 µl	50 µl
10 kb primers	4 µM each	50 µl	50 µl
DNA size standard		200 µl	400 µl
DMSO		0.5 ml	0.5 ml

## 3. Guidelines for Phusion™ DNA Polymerase

Phusion DNA Polymerase (2U/µl) is provided with 5x Phusion HF Buffer and 5x Phusion GC Buffer. Both buffers contain 1.5 mM MgCl<sub>2</sub> at final reaction concentrations. Separate tubes of DMSO and 50 mM MgCl<sub>2</sub> solutions are provided for further optimization.

### 3.1 Basic reaction conditions for DNA amplification

Carefully mix and centrifuge all tubes before opening to improve recovery. PCR reactions should be set up on ice. Phusion DNA Polymerase should be pipetted carefully and gently as the high glycerol content (50 %) in the storage buffer may otherwise lead to pipetting errors. It is critical that the Phusion DNA Polymerase is the last component added to the PCR mixture, since the enzyme exhibits 3'→5' exonuclease activity that can degrade primers in the absence of dNTPs.

**Table 1. Pipetting instructions (add items in this order).**

Component	Volume / 50 µl reaction	Volume / 20 µl reaction	Final conc.
H <sub>2</sub> O	add to 50 µl	add to 20 µl	
5x Phusion HF Buffer*	10 µl	4 µl	1x
10 mM dNTPs	1 µl	0.4 µl	200 µM each
primer A**	x µl	x µl	0.5 µM
primer B**	x µl	x µl	0.5 µM
template DNA	x µl	x µl	
(DMSO***, optional)	(1.5 µl)	(0.6 µl)	(3 %)
Phusion DNA Polymerase	0.5 µl	0.2 µl	0.02U/µl

\* Optionally 5x Phusion GC Buffer can be used, see section 4.2 for details.

\*\* The recommendation for final primer concentration is 0.5 µM, but it can be varied in a range of 0.2-1.0 µM, if needed.

\*\*\* Addition of DMSO is recommended for GC-rich amplicons. DMSO is not recommended for amplicons with very low GC % or amplicons that are >20 kb.

### 3.2 Cycling conditions

Due to the novel nature of Phusion DNA Polymerase, optimal reaction conditions may differ from standard enzyme protocols. Phusion DNA Polymerase tends to work better at elevated denaturation and annealing temperatures due to higher salt concentrations in its buffer. Following the guidelines will ensure optimal enzyme performance.

**Table 2. Cycling instructions.**

Cycle step	2-step protocol		3-step protocol		Cycles
	Temp.	Time	Temp.	Time	
Initial denaturation	98°C	30 s	98°C	30 s	1
Denaturation	98°C	5-10 s	98°C	5-10 s	25-35
Annealing (see 5.3)	-	-	X°C	10-30 s	
Extension (see 5.4)	72°C	15-30 s/1 kb	72°C	15-30 s/1 kb	
Final extension	72°C 4°C	5-10 min hold	72°C 4°C	5-10 min hold	1

## 4. Notes about Reaction Components

### 4.1 Enzyme

The optimal amount of enzyme depends on the amount of template and the length of the PCR product. Usually 1 unit of Phusion DNA Polymerase per 50 µl reaction volume gives good results, but optimal amounts could range from 0.5-2 units per 50 µl reaction depending on amplicon length and difficulty. **Do not exceed 2 U/50 µl (0.04 U/µl), especially for amplicons that are > 5kb.**

When cloning fragments amplified with Phusion DNA Polymerase, blunt end cloning is recommended. If TA cloning is required, it can be performed by adding A overhangs to the blunt PCR product with e.g. DyNAzyme™ II DNA Polymerase (F-501). However, before adding the overhangs it is very important to remove all the Phusion DNA Polymerase by purifying the PCR product carefully. Any remaining Phusion DNA Polymerase will degrade the A overhangs, thus creating blunt ends again. A detailed protocol for TA cloning of PCR products amplified with any of the Phusion DNA polymerases can be found on Finnzymes website ([www.finnzymes.com](http://www.finnzymes.com)).

### 4.2 Buffers

Two buffers are provided with the enzyme: 5x Phusion HF Buffer and 5x Phusion GC Buffer. The error rate of Phusion DNA Polymerase in HF Buffer ( $4.4 \times 10^{-7}$ ) is lower than that in GC Buffer ( $9.5 \times 10^{-7}$ ). Therefore, the HF Buffer should be used as the default buffer for high-fidelity amplification. However, GC Buffer can improve the performance of Phusion DNA Polymerase on some difficult or long templates, i.e. GC-rich templates or those with complex secondary structures.

### 4.3 Mg<sup>2+</sup> concentration and dNTP concentration

Concentration of Mg<sup>2+</sup> is critical since Phusion DNA Polymerase is a magnesium dependent enzyme. However, excessive Mg<sup>2+</sup> stabilizes the DNA double strand and prevents complete denaturation of DNA. Excess Mg<sup>2+</sup> can also stabilize spurious annealing of primers to incorrect template sites and decrease specificity. Conversely, inadequate Mg<sup>2+</sup> could lead to lower product yield. The optimal Mg<sup>2+</sup> concentration will also depend on the dNTP concentration, the specific template DNA and the sample buffer composition. In general, the optimal Mg<sup>2+</sup> concentration is 0.5 - 1 mM over the total dNTP concentration for standard PCR. If the primers and/or template contain chelators such as EDTA or EGTA, the apparent Mg<sup>2+</sup>-optimum may be shifted to higher concentrations. If further optimization is needed, increase Mg<sup>2+</sup>-concentration in 0.2 mM steps.

High quality dNTPs (e.g. F-560) should be used for optimal performance with Phusion DNA Polymerase. Use of dUTP, dITP, or dUTP-derivatives or analogues is not recommended. Due to the increased processivity of Phusion DNA Polymerase there is no advantage in increasing dNTP concentrations. For optimal results always use 200 µM of each dNTP.

### 4.4 Template

General guidelines are: 1 pg - 10 ng / 50 µl reaction with low complexity DNA (e.g. plasmid, lambda or BAC DNA); 50-250 ng / 50 µl reaction with high complexity genomic DNA. If cDNA synthesis reaction mixture is used as a source of template, the volume used should not exceed 10 % of the final PCR reaction volume.

### 4.5 PCR additives

The recommended reaction conditions for GC-rich templates include 3 % DMSO as a PCR additive, which aids in the denaturing of templates with high GC contents. For further optimization DMSO should be varied in 2 % increments. In some cases DMSO may also be required for supercoiled plasmids to relax for denaturation. Other PCR additives such as formamide, glycerol, and betaine are also compatible with Phusion DNA Polymerase.

If high DMSO concentration is used, the annealing temperature must be lowered, as DMSO decreases the melting point of the primers. It has been reported that 10 % DMSO decreases the annealing temperature by 5.5-6.0°C.<sup>2</sup>

## 5. Notes about Cycling Conditions

### 5.1 Initial denaturation

Denaturation should be done at 98°C. Due to the high thermostability of Phusion DNA Polymerase even higher denaturation temperatures can be used. We recommend 30 seconds initial denaturation at 98°C for most templates. Some templates may require longer initial denaturation and the length of the initial denaturation time can be extended up to 3 minutes.

### 5.2 Denaturation

Keep the denaturation as short as possible. Usually 5-10 seconds at 98°C is enough for most templates. **Note:** The denaturation time and temperature may vary depending on the ramp rate and temperature control mode of the cycler.

### 5.3 Primer annealing

The Phusion DNA Polymerase has the ability to stabilize primer-template hybridization. As a basic rule, for primers > 20nt, anneal for 10 - 30 seconds at T<sub>m</sub> +3°C of the lower T<sub>m</sub> primer. The T<sub>m</sub>'s should be calculated with the nearest-neighbor method<sup>3</sup> as results from primer T<sub>m</sub> calculations can vary significantly depending on the method used. For primers ≤ 20nt, use an annealing temperature equal to the T<sub>m</sub> of the lower T<sub>m</sub> primer. If necessary, use a temperature gradient to find the optimal annealing temperature for each template-primer pair combination. The annealing gradient should extend up to the extension temperature (two-step PCR). Two-step cycling without annealing step is also recommended for high T<sub>m</sub> primer pairs. Instructions for T<sub>m</sub> calculation and a link to a calculator using the nearest-neighbor method can be found on Finnzymes' website ([www.finnzymes.com](http://www.finnzymes.com)).

### 5.4 Extension

The extension should be performed at 72°C. Extension time depends on amplicon length and complexity. For low complexity DNA (e.g. plasmid, lambda or BAC DNA) use extension time 15 seconds per 1kb. For high complexity genomic DNA 30 seconds per 1kb is recommended. For some cDNA templates, the extension time can be increased up to 40 seconds per 1 kb to obtain optimal results.

## 6. Amplifying Control Template

### 6.1 Reaction conditions

Table 3. Pipetting instructions for control reactions.

Component	Volume / 50 µl reaction	Volume / 20 µl reaction	Final conc.
H <sub>2</sub> O	34 µl	13.6 µl	
5x Phusion HF Buffer	10 µl	4 µl	1x
10 mM dNTPs	1 µl	0.4 µl	200 µM each
Primers*	2.5 µl	1 µl	0.2 µM
Control template DNA	2 µl	0.8 µl	
Phusion DNA Polymerase	0.5 µl	0.2 µl**	0.02U/µl

\* Either the 1.3 kb primer set or 10 kb primer set.

\*\* Dilution of polymerase should be made to 1x reaction buffer to avoid pipetting errors.

### 6.2 Cycling conditions

A separate cycling protocol is given for both 1.3 kb and 10 kb control amplicons. Alternatively, both control reactions can be amplified simultaneously using the 10 kb cycling protocol.

Table 4. Cycling conditions for 1.3 kb control fragment (2-step protocol).


Cycle step	Temp.	Time	Number of cycles
Initial denaturation	98°C	1 min	1
Denaturation	98°C	5 s	25
Annealing / Extension	72°C	20 s	
Final extension	72°C	10 min	1
	10°C	hold	

**Table 5. Cycling conditions for 10 kb fragment (3-step protocol). This program can also be used if both control reactions are amplified simultaneously.**

Cycle step	Temp.	Time	Number of cycles
Initial denaturation	98°C	1 min	1
Denaturation	98°C	5 s	25
Annealing	60°C	15 s	
Extension	72°C	2 min 30 s	
Final extension	72°C 10°C	10 min hold	1

The cycling protocols above are recommendations. If you wish to run these controls together or with your experimental samples, please note that the controls have been shown to work in a variety of conditions. The 1.3 kb control has been successfully amplified with both 2- and 3-step protocols with extension times ranging from 15 s to 5 min, and cycle numbers ranging from 20 to 30. The 10 kb control has been successfully amplified with 3-step protocol with extension times ranging from 2 min to 5 min, and cycle numbers ranging from 20 to 30.

### 6.3 Analysis of the control reactions



In the image on the left both control reactions have been run on an ethidium bromide stained agarose gel (1% SeaKem LE agarose in TAE buffer). For this run 15 µl of loading dye was added to the 50 µl control PCR reactions, and 5 µl of the resulting mixtures were loaded on the gel.

Lane 1. DNA size standard  
Lane 2. 1.3 kb control amplicon  
Lane 3. 10 kb control amplicon  
Lane 4. DNA size standard

After running your control reactions on a gel, compare the results to the image on the left to check for specificity and efficiency of the reactions.

## 7. Troubleshooting

No product at all or low yield
<ul style="list-style-type: none"> <li>Repeat and make sure that there are no pipetting errors.</li> <li>Use fresh high quality dNTPs. Do not use dNTP mix that contains dUTP.</li> <li>Titrate template amount.</li> <li>Template DNA may be damaged. Use carefully purified template.</li> <li>Increase extension time.</li> <li>Increase the number of cycles.</li> <li>Optimize annealing temperature.</li> <li>Optimize enzyme concentration.</li> <li>Titrate DMSO (2-8 %) in the reaction (see section 4.5).</li> <li>Denaturation temperature may be too low. Optimal denaturation temperature for most templates is 98°C or higher.</li> <li>Optimize denaturation time.</li> <li>Check the purity and concentration of the primers.</li> <li>Check primer design.</li> <li>Try using the alternative GC buffer (see section 4.2).</li> <li>If DNA is not carefully purified, inhibitors may be present - decrease the amount of DNA.</li> </ul>
Non-specific products - High molecular weight smears
<ul style="list-style-type: none"> <li>Reduce enzyme concentration (see section 4.1).</li> <li>Shorten extension time (see section 5.4).</li> <li>Titrate template amount.</li> <li>Reduce the total number of cycles.</li> <li>Increase annealing temperature or try 2-step protocol (see section 5.3).</li> <li>Vary denaturation temperature (see section 5.2).</li> <li>Optimize Mg<sup>2+</sup>-concentration.</li> <li>Decrease primer concentration.</li> </ul>
Non-specific products - Low molecular weight discrete bands
<ul style="list-style-type: none"> <li>Increase annealing temperature (see section 5.3).</li> <li>Shorten extension time (see section 5.4).</li> <li>Decrease enzyme concentration.</li> <li>Optimize Mg<sup>2+</sup>-concentration.</li> <li>Titrate template amount.</li> <li>Decrease primer concentration.</li> <li>Design new primers.</li> </ul>

## 8. Component Specifications

### 8.1 Phusion™ High-Fidelity DNA Polymerase (F-530)

Thermostable Phusion DNA Polymerase is purified from an *E.coli* strain expressing the cloned Phusion DNA Polymerase gene. Phusion DNA Polymerase possesses the following activities: 5'→3' DNA polymerase activity and 3'→5' exonuclease activity. Phusion DNA Polymerase is purified free of contaminating endo- and exonucleases.

**Storage buffer:** 20 mM Tris-HCl (pH 7.4 at 25°C), 0.1 mM EDTA, 1 mM DTT, 100 mM KCl, stabilizers, 200 µg/ml BSA and 50 % glycerol.

**Unit definition:** One unit is defined as the amount of enzyme that will incorporate 10 nmoles of dNTPs into acid-insoluble form at 74°C in 30 minutes under the stated assay conditions.

**Unit assay conditions:** Incubation buffer: 25 mM TAPS-HCl, pH 9.3 (at 25°C), 50 mM KCl, 2 mM MgCl<sub>2</sub>, 1 mM β-mercaptoethanol, 100 µM dCTP, 200 µM each dATP, dGTP, dTTP. Incubation procedure: 20 µg activated calf thymus DNA and 0.5 µCi [α-<sup>32</sup>P] dCTP are incubated with 0.1 units of DNA polymerase in 50 µl incubation buffer at 74°C for 10 minutes. The amount of incorporated dNTPs is determined by trichloroacetic acid precipitation.

**DNA amplification assay:** Performance in PCR is tested by the amplification of 7.5 kb genomic DNA and 20 kb lambda DNA.

**Exonuclease activity:** Incubation of 10 U for 4 hours at 72°C in 50 µl assay buffer with 1 µg sonicated <sup>3</sup>H ssDNA (2x10<sup>5</sup> cpm/µg) released < 1 % of radioactivity.

**Endonuclease assay:** No endonuclease activity is observed after incubation of 10 U of DNA polymerase with 1 µg of λ DNA in assay buffer at 72°C for 4 hours.

### 8.2 5x Phusion™ HF Buffer (F-518)

The 5x Phusion HF Buffer contains 7.5 mM MgCl<sub>2</sub>, which provides 1.5 mM MgCl<sub>2</sub> in final reaction conditions.

### 8.3 5x Phusion™ GC Buffer (F-519)

The 5x Phusion GC Buffer contains 7.5 mM MgCl<sub>2</sub>, which provides 1.5 mM MgCl<sub>2</sub> in final reaction conditions.

### 8.4 dNTP mix (F-560)

The dNTP mix is a premixed ready-to-use solution consisting of the following compounds: dATP, dGTP, dCTP and dTTP dissolved in H<sub>2</sub>O at 10 mM each.

### 8.5 50 mM MgCl<sub>2</sub> solution (F-510MG)

Both Phusion Buffers supply 1.5 mM MgCl<sub>2</sub> at final reaction conditions. If higher MgCl<sub>2</sub> concentrations are desired, use this 50 mM MgCl<sub>2</sub> solution to increase the MgCl<sub>2</sub> titer. Using the following equation you can calculate the volume of 50 mM MgCl<sub>2</sub> needed to attain the final MgCl<sub>2</sub> concentration:

[desired mM Mg] - [1.5 mM] = µl to add to a 50 µl reaction

Because the PCR reactions can be quite sensitive to changes in the MgCl<sub>2</sub> concentration, it is recommended that the 50 mM MgCl<sub>2</sub> stock solution is diluted 1:5 (to 10 mM) to minimize pipetting errors.

### 8.6 Lambda DNA control template (F-304K)

The control template is bacteriophage lambda DNA (GenBank access number J02459, 48 502 bp). The concentration is 0.5 ng/µl in TE buffer.

### 8.7 1.3 kb control primer mix (F-535)

This component is a mix of primers in H<sub>2</sub>O for a 1.3 kb fragment of lambda DNA. Each primer concentration is 4 µM. Below are the details for the primers:

Primer #1 (27-mer)  
5'-GTC ACC AGT GCA GTG CTT GAT AAC AGG-3'  
Melting point: 71.0°C  
Coordinates in lambda DNA: 30 006 - 30 032

Primer #2 (28-mer)  
5'-GAT GAC GCA TCC TCA CGA TAA TAT CCG G-3'  
Melting point: 73.2°C  
Coordinates in lambda DNA: 31 325 - 31 352

### 8.8 10 kb control primer mix (F-536)

This component is a mix of primers in H<sub>2</sub>O for a 10 kb fragment of lambda DNA. Each primer concentration is 4 µM.


Primer #1 (22-mer)  
5'-CAG TGC AGT GCT TGA TAA CAG G-3'  
Melting point: 63.5°C  
Coordinates in lambda DNA: 30 011 - 30 032

Primer #2 (20-mer)  
5'-GTA GTG CGC GTT TGA TTT CC-3'  
Melting point: 63.3°C  
Coordinates in lambda DNA: 40 024 - 40 043

### 8.9 Ready-to-use DNA size and mass standard (F-303SD)

This size and mass standard is a mix of lambda DNA *Hind*III digest and bacteriophage φX174 DNA *Hae*III digest, each at 50 ng/µl (100 ng/µl total). It is supplied in 8 mM Tris-HCl (pH 8.0), 12 mM EDTA, 12 % glycerol and 0.012 % (w/v) bromophenol blue dye.

The DNA standard solution contains 19 fragments of the following sizes and mass amounts (per 10 µl):



Fragment	Base pairs	DNA amount ng/10 µl
1	23 130	238
2	9 416	97
3	6 557	68
4	4 361	45
5	2 322	24
6	2 027	21
7	1 353	126
8	1 078	100
9	872	81
10	603	56
11	564*	6
12	310	29
13a	281	26
13b	271	25
14	234	22
15	194	18
16a	125*	1
16b	118	11
17	72	7

**Note:** The cohesive areas of fragments 1 and 4 can be separated by heating at 65°C for 5 minutes. For daily use the marker can be stored at +4°C (at least one month). The marker is stable at -20°C for at least one year.

\* Due to the low amount of DNA these bands are almost invisible.

### 8.10 Dimethyl sulfoxide DMSO, 100 % (F-515)

**Note:** The freezing point of DMSO is 18-19°C, so it does not melt on ice.

## 9. References

1. Frey & Suppmann, (1995) *Biochemica* 34-35.
2. Chester & Marshak, (1993) *Analytical Biochemistry* 209, 284-290.
3. Breslauer *et al.*, (1986) *PNAS* 83, 3746-3750.

### Storage and shipping

Phusion High-Fidelity PCR Kit is shipped on gel ice. Upon arrival, store the components at -20°C. The Phusion High-Fidelity PCR Kit is stable for one year from the packaging date when stored and handled properly.

### Warranty

Finnzymes Oy warrants that its products will meet the specifications stated on the technical data section of the data sheets, and Finnzymes Oy agrees to replace the products free of charge if the products do not conform to the specifications. Notice for replacement must be given within 60 days of receipt. In consideration of the above commitments by Finnzymes Oy, the buyer agrees to and accepts the following conditions:

- That this warranty is in lieu of all other warranties, express or implied;
- That **ALL WARRANTIES OF MERCHANT ABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE ARE HEREBY EXCLUDED AND WAIVED**;
- That the buyer's sole remedy shall be to obtain replacement of the product free of charge from Finnzymes Oy; and
- That this remedy is in lieu of all other remedies or claims for damages, consequential or otherwise, which the buyer may have against Finnzymes Oy.

### Exclusive terms of sale

Finnzymes Oy does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorised officer of Finnzymes Oy. Prices are subject to change without notice.

### Recommended guidelines for safe use of the products

Finnzymes Oy recommends that the buyer and other persons using the products follow the N.I.H. guidelines published in the Federal Register, Volume 41, No. 131, July 7, 1976, and any amendments thereto. Finnzymes Oy disclaims any and all responsibility for any injury or damage which may be caused by the failure of the buyer or any other person to follow said guidelines.

### Research use only

Since these products are intended for research purposes by qualified persons, the Environmental Protection Agency does not require us to supply Premanufacturing Notice.

### Notice to user

The information presented here is accurate and reliable to the best of our knowledge and belief, but is not guaranteed to be so. Nothing herein is to be construed as recommending any practice or any product in violation of any patent or in violation of any law or regulation. It is the user's responsibility to determine for himself or herself the suitability of any material and/or procedure for a specific purpose and to adopt such safety precautions as may be necessary.

### Phusion™ DNA Polymerases Notice to Purchaser

Limited license (proofreading DNA polymerases). The purchase price of this product includes a limited, non-transferable license under U.S. and foreign patents (5,500,363 and 5,352,778) owned by New England Biolabs, Inc. to use this product. No other license under these patents is conveyed expressly or by implication to the purchaser by the purchase of this product. The purchase price of this product includes a limited, non-transferable license under U.S. and foreign patents owned by BIO-RAD Laboratories, Inc., to use this product. No other license under these patents is conveyed expressly or by implication to the purchaser by the purchase of this product.

The quality system of Finnzymes Oy is certified according to standard SFS-EN ISO9001:2000.

Phusion is a trademark of Finnzymes Oy.

Version 1.5, December 2007