

Phusion® Hot Start High-Fidelity DNA Polymerase

Product codes: F-540S, 100 U
F-540L, 500 U

Stable for one year from the assay date. Store at –20°C.

1. Introduction

Finnzymes' Phusion® Hot Start High-Fidelity DNA Polymerase offers superior performance for all PCR applications. A unique processivity-enhancing domain makes this *Pyrococcus*-like proofreading enzyme extremely processive, accurate and rapid. The error rate of Phusion Hot Start DNA Polymerase is equal to that of Phusion DNA Polymerase (4.4×10^{-7} in Phusion HF-buffer) when determined with a modified *Ia*cl-based method¹. It is approximately 50-fold lower than that of *Thermus aquaticus* DNA polymerase and 6-fold lower than that of *Pyrococcus furiosus* DNA polymerase. Phusion Hot Start High-Fidelity DNA Polymerase is capable of amplifying long amplicons such as the 7.5 kb genomic and 20 kb λ DNA used in Finnzymes' quality control assays.

Phusion Hot Start DNA Polymerase combines the DNA polymerase and a reversibly bound, specific Affibody® protein^{2,3}, which inhibits the DNA polymerase activity at ambient temperatures, thus preventing the amplification of non-specific products. In addition, the Affibody ligand inhibits the 3'→5' exonuclease activity of the polymerase, preventing degradation of primers and template DNA during reaction setup. At polymerization temperatures, the Affibody molecule is released, rendering the polymerase fully active. Phusion Hot Start DNA Polymerase does not require any separate activation step in the PCR protocol.

Phusion Hot Start DNA Polymerase possesses the following activities: 5'→3' DNA polymerase activity and 3'→5' exonuclease activity. It generates blunt ends in the amplification products.

Important notes

- Use Phusion DNA Polymerase at 0.5–1.0 U per 50 μ l reaction volume. Do not exceed 2 U/50 μ l. (See 4.1)
- Use 15–30 s/kb for extension. Do not exceed 1 min/kb. (See 5.4)
- Use 98°C for denaturation. (See 5.1 & 5.2)
- Anneal at $T_m+3^\circ\text{C}$ (> 20nt) or use 2-step protocol. (See 5.3)
- Use 200 μ M of each dNTP. Do not use dUTP. (See 4.3)
- Note: Phusion DNA Polymerases produce blunt end DNA products.

2. Package information

F-540S	100 U (2 U/μl) Material provided: Phusion® Hot Start DNA Polymerase 100 U (50 μ l), 5x Phusion® HF Buffer (2 x 1.5 ml), 5x Phusion® GC Buffer (1.5 ml), DMSO (500 μ l) and 50 mM MgCl ₂ solution (1.5 ml).
F-540L	500 U (2 U/μl) Material provided: Phusion® Hot Start DNA Polymerase 500 U (250 μ l), 5x Phusion® HF Buffer (6 x 1.5 ml), 5x Phusion® GC Buffer (2 x 1.5 ml), DMSO (500 μ l) and 50 mM MgCl ₂ solution (2 x 1.5 ml).

Reaction buffer: 5x Phusion HF Buffer and 5x Phusion GC Buffer both contain 7.5 mM MgCl₂.

Material safety data sheet (MSDS) is available at www.finnzymes.com.

3. Guidelines for using Phusion® Hot Start DNA Polymerase

Carefully mix and centrifuge all tubes before opening to ensure homogeneity and improve recovery. When using Phusion Hot Start DNA Polymerase, it is not necessary to perform the PCR setup on ice. Prepare a master mix for the appropriate number of samples to be amplified. The DNA polymerase should be pipetted carefully and gently as the high glycerol content (50 %) in the storage buffer may otherwise lead to pipetting errors.

Protocols optimized for Phusion® DNA Polymerase can be applied to Phusion Hot Start DNA Polymerase reactions. However, use an annealing temperature of 60°C or higher (see section 5.3). Due to the novel nature of Phusion Hot Start DNA Polymerase, the optimal reaction conditions may differ from PCR protocols for standard DNA polymerases. Due to the high salt concentration in the reaction buffer, Phusion Hot Start DNA Polymerase tends to work better at elevated denaturation and annealing temperatures. Please pay special attention to the conditions listed below when running your reactions. Following the guidelines will ensure optimal enzyme performance.

Table 1. Pipetting instructions (add items in this order).

Component	50 μ l reaction	20 μ l reaction	Final conc.
H ₂ O	add to 50 μ l	add to 20 μ l	
5x Phusion® HF Buffer*	10 μ l	4 μ l	1x
10 mM dNTPs	1 μ l	0.4 μ l	200 μ M each
primer A**	x μ l	x μ l	0.5 μ M
primer B**	x μ l	x μ l	0.5 μ M
template DNA	x μ l	x μ l	
(DMSO***, optional)	(1.5 μ l)	(0.6 μ l)	(3 %)
Phusion® Hot Start DNA Polymerase (2 U/ μ l)	0.5 μ l	0.2 μ l	0.02 U/ μ l

* Optionally 5x Phusion GC Buffer can be used. See section 4.2. for details.

** The recommendation for final primer concentration is 0.5 μ M, but it can be varied in a range of 0.2–1.0 μ M if needed.

*** Addition of DMSO is recommended for GC-rich amplicons. DMSO is not recommended for amplicons with very low GC % or amplicons that are >20 kb.

Table 2. Cycling instructions.

Cycle step	2-step protocol		3-step protocol		Cycles
	Temp.	Time	Temp.	Time	
Initial denaturation	98°C	30 s	98°C	30 s	1
Denaturation	98°C	5–10 s	98°C	5–10 s	25–35
Annealing (see 5.3)	–	–	X°C	10–30 s	
Extension (see 5.4)	72°C	15–30 s/kb	72°C	15–30 s/kb	
Final extension	72°C 4°C	5–10 min hold	72°C 4°C	5–10 min hold	1

4. Notes about reaction components

4.1 Enzyme

The optimal amount of enzyme depends on the amount of template and the length of the PCR product. Usually 1 unit of Phusion Hot Start DNA Polymerase per 50 µl reaction volume gives good results, but the optimal amount can range from 0.5 to 2 units per 50 µl reaction depending on the amplicon length and difficulty. **Do not exceed 2 U/50 µl (0.04 U/µl), especially for amplicons that are > 5 kb.**

When cloning fragments amplified with Phusion Hot Start DNA Polymerase, blunt end cloning is recommended. If TA cloning is required, it can be performed by adding A overhangs to the blunt PCR product with DyNAzyme™ II DNA Polymerase (F-501), for example. However, before adding the overhangs it is very important to remove all Phusion Hot Start DNA Polymerase by purifying the PCR product carefully. Any remaining Phusion Hot Start DNA Polymerase will degrade the A overhangs, creating blunt ends again. A detailed protocol for TA cloning of PCR fragments amplified with any of the Phusion DNA Polymerases can be found on Finnzymes' website (www.finnzymes.com).

4.2 Buffers

Two buffers are provided with the enzyme: 5x Phusion HF Buffer (F-518) and 5x Phusion GC Buffer (F-519). The error rate of Phusion Hot Start DNA Polymerase in HF Buffer (4.4×10^{-7}) is lower than that in GC Buffer (9.5×10^{-7}). Therefore, HF Buffer should be used as the default buffer for high-fidelity amplification. However, GC Buffer can improve the performance of Phusion Hot Start DNA Polymerase on some difficult or long templates, such as GC-rich templates or those with complex secondary structures. For applications such as microarray or DHPLC, where the DNA templates need to be free of detergents, detergent-free reaction buffers (F-520, F-521) are available for Phusion DNA Polymerases.

4.3 Mg²⁺ and dNTP

The concentration of Mg²⁺ is critical since Phusion Hot Start DNA Polymerase is a magnesium-dependent enzyme. Excessive Mg²⁺ stabilizes the DNA double strand and prevents complete denaturation of DNA. Excess Mg²⁺ can also stabilize spurious annealing of primers to incorrect template sites and decrease specificity. Conversely, inadequate Mg²⁺ may lead to lower product yield. The optimal Mg²⁺ concentration also depends on the dNTP concentration, the specific template DNA and the sample buffer composition. In general, the optimal Mg²⁺ concentration is 0.5 to 1 mM over the total dNTP concentration for standard PCR. If the primers and/or template

contain chelators such as EDTA or EGTA, the apparent Mg²⁺ optimum may be shifted to higher concentrations. If further optimization is needed, increase Mg²⁺ concentration in 0.2 mM steps.

High quality dNTPs (e.g. F-560S/L) should be used for optimal performance with Phusion Hot Start DNA Polymerase. The polymerase cannot read dUTP-derivatives or dITP in the template strand so the use of these analogues or primers containing them is not recommended. Due to the high processivity of Phusion Hot Start DNA Polymerase there is no advantage of increasing dNTP concentrations. For optimal results always use 200 µM of each dNTP.

4.4 Template

General guidelines for low complexity DNA (e.g. plasmid, lambda or BAC DNA) are: 1 pg–10 ng per 50 µl reaction volume. For high complexity genomic DNA, the amount of DNA template should be 50–250 ng per 50 µl reaction volume. If cDNA synthesis reaction mixture is used as a source of template, the volume of the template should not exceed 10 % of the final PCR reaction volume.

4.5 PCR additives

The recommended reaction conditions for GC-rich templates include 3 % DMSO as a PCR additive, which aids in the denaturing of templates with high GC content. For further optimization DMSO should be increased in 2 % steps. In some cases DMSO may also be required for supercoiled plasmids to relax for denaturation. Other PCR additives such as formamide (up to 3 %), glycerol and betaine are also compatible with Phusion Hot Start DNA Polymerase.

If high DMSO concentration is used, the annealing temperature must be decreased, as DMSO affects the melting point of the primers. It has been reported that 10 % DMSO decreases the annealing temperature by 5.5–6.0°C⁴.

5. Notes about cycling conditions

5.1 Initial denaturation

Denaturation should be performed at 98°C. Due to the high thermostability of Phusion Hot Start DNA Polymerase even higher than 98°C denaturation temperatures can be used. We recommend a 30-second initial denaturation at 98°C for most templates. Some templates may require longer initial denaturation, and the length of the initial denaturation time can be extended up to 3 minutes.

5.2 Denaturation

Keep the denaturation as short as possible. Usually 5–10 seconds at 98°C is enough for most templates. **Note:** the denaturation time and temperature may vary depending on the ramp rate and temperature control mode of the cyclor.

5.3 Primer annealing

With Phusion Hot Start DNA Polymerase, use primers with T_m 60°C or higher. Typically the length of such primers is 20 nt or more. The T_m values should be calculated with the nearest-neighbor method⁵ as results from primer T_m calculations can vary significantly depending on the method used. Instructions

for T_m calculation and a link to a calculator using the nearest-neighbor method can be found on Finnzymes' website (www.finnzymes.com). As a basic rule, for primers >20 nt, anneal for 10–30 seconds at a T_m +3°C of the lower T_m primer. For primers ≤ 20 nt, use an annealing temperature equal to the T_m of the lower T_m primer. If necessary, use a temperature gradient to find the optimal annealing temperature for each template-primer pair combination. The annealing gradient should extend up to the extension temperature (two-step PCR). A 2-step protocol is recommended when primer T_m values are at least 69°C (> 20 nt) or 72°C (≤ 20 nt) when calculated with Finnzymes' T_m calculator. In the 2-step protocol the combined annealing/extension step should be performed at 72°C even when the primer T_m is > 72°C.

5.4 Extension

The extension should be performed at 72°C. The extension time depends on the length and complexity of the amplicon. For low complexity DNA (e.g. plasmid, lambda or BAC DNA) use an extension time of 15 seconds per 1 kb. For high complexity genomic DNA, 30 seconds per 1 kb is recommended. For some cDNA templates, the extension time can be increased up to 40 seconds per 1 kb to obtain optimal results.

6. Troubleshooting

No product at all or low yield

- Repeat the PCR and make sure that there are no pipetting errors.
- Use primers with T_m 60°C or higher. T_m values can be increased by lengthening the primers. Preferable primer length is 20 nt or longer (see section 5.3).
- Use Finnzymes' T_m calculator (www.finnzymes.com).
- Use fresh high-quality dNTPs.
- Do not use dNTP mix or primers that contain dUTP or dITP.
- Sample concentration may be too low. Use more template.
- Template DNA may be damaged. Use carefully purified template.
- Increase extension time.
- Increase the number of cycles.
- Decrease annealing temperature.
- Optimize enzyme concentration.
- Titrate DMSO (2–8 %) in the reaction (see section 4.5).
- Denaturation temperature may be too low. Optimal denaturation temperature for most templates is 98°C or higher.
- Denaturation time may be too long or too short. Optimize denaturation time.
- Check the purity and concentration of the primers.
- Check primer design.
- Try using the alternative GC Buffer (see section 4.2).

Non-specific products - High molecular weight smears

- decrease enzyme concentration (see section 4.1).
- Decrease extension time (see section 5.4).
- Reduce the total number of cycles.
- Increase annealing temperature or try 2-step protocol (see section 5.3).
- Vary denaturation temperature (see section 5.2).
- Optimize Mg²⁺ concentration (see section 4.3).
- Reduce primer concentration.

Non-specific products - Low molecular weight discrete bands

- Increase annealing temperature (see section 5.3).
- Shorten extension time (see section 5.4).
- Reduce enzyme concentration (see section 4.1).
- Optimize Mg²⁺ concentration (see section 4.3).
- Titrate template amount.
- Decrease primer concentration.
- Design new primers.

7. Component specifications

7.1 Phusion® Hot Start High-Fidelity DNA Polymerase (F-540)

Thermostable Phusion DNA Polymerase is isolated and purified from an *E. coli* strain expressing the cloned Phusion DNA Polymerase gene. Phusion DNA Polymerase possesses the following activities: 5'→3' DNA polymerase activity and 3'→5' exonuclease activity. The Affibody ligand is isolated and purified from an *E. coli* strain expressing the cloned Affibody-encoding gene. Phusion Hot Start DNA Polymerase is free of contaminating endo- and exonucleases.

Storage buffer: 20 mM Tris-HCl (pH 7.4 at 25°C), 0.1 mM EDTA, 1 mM DTT, 100 mM KCl, stabilizers, 200 µg/ml BSA and 50 % glycerol.

Unit definition: One unit is defined as the amount of enzyme that will incorporate 10 nmoles of dNTPs into acid-insoluble form at 74°C in 30 minutes under the stated assay conditions.

Unit assay conditions: Incubation buffer: 25 mM TAPS-HCl, pH 9.3 (at 25°C), 50 mM KCl, 2 mM MgCl₂, 1 mM β-mercaptoethanol, 100 µM dCTP, 200 µM each dATP, dGTP, dTTP.

Incubation procedure: 20 µg activated calf thymus DNA and 0.5 µCi [α -³²P] dCTP are incubated with 0.1 units of DNA polymerase in 50 µl incubation buffer at 74°C for 10 minutes. The amount of incorporated dNTPs is determined by trichloroacetic acid precipitation.

DNA amplification test: Performance in PCR is tested by the amplification of 2.3 and 7.5 kb fragments of genomic DNA and 20 kb fragment of λ DNA.

Exonuclease contamination assay: Incubation of 10 U for 4 hours at 72°C in 50 µl assay buffer with 1 µg sonicated [³H] DNA (2 × 10⁵ cpm/µg) released < 1 % of radioactivity.

Endonuclease contamination assay: No endonuclease activity was observed after incubation of 10 U of DNA polymerase with 1 µg of λ DNA in assay buffer at 72°C for 4 hours.

7.2 5x Phusion® HF Buffer (F-518)

The 5x Phusion HF Buffer contains 7.5 mM MgCl₂, which provides 1.5 mM MgCl₂ in final reaction conditions.

7.3 5x Phusion® GC Buffer (F-519)

The 5x Phusion GC Buffer contains 7.5 mM MgCl₂, which provides 1.5 mM MgCl₂ in final reaction conditions.

Caution: Repeated freezing and thawing of the buffer can result in the precipitation or accumulation of MgCl₂ in insoluble form. For consistent results, heat the buffer to 90°C for 10 min and vortex prior to use if needed, or store refrigerated.

7.4 50 mM MgCl₂ Solution (F-510MG)

Both Phusion Buffers supply 1.5 mM MgCl₂ at final reaction conditions. If higher MgCl₂ concentrations are desired, use a 50 mM MgCl₂ solution to increase the MgCl₂ titer. Using the following equation, you can calculate the volume of 50 mM MgCl₂ needed to attain the final MgCl₂ concentration: [desired mM Mg] – [1.5 mM] = μ l to add to a 50 μ l reaction.

For example to increase the MgCl₂ concentration to 2.0 mM, add 0.5 μ l of the 50 mM MgCl₂ solution. Because the PCR reactions can be quite sensitive to changes in the MgCl₂ concentration, it is recommended that the 50 mM MgCl₂ stock solution is diluted 1:5 (to 10 mM) to minimize pipetting errors.

8. References

1. Frey M. & Suppmann B. (1995) *Biochemica* 2: 34–35.
2. Nord K. *et al.* (1997) *Nature Biotechnol.* 15: 772–777.
3. Wikman M. *et al.* (2004) *Protein Eng. Des. Sel.* 17: 455–462.
4. Chester N. & Marshak D.R. (1993) *Anal. Biochem.* 209: 284–290.
5. Breslauer K.J. *et al.* (1986) *PNAS* 83: 3746–3750.

Shipping and storage

Phusion Hot Start DNA Polymerase is shipped on gel ice. Upon arrival, store the components at –20°C. The Phusion Hot Start DNA Polymerase is stable for one year from the assay date when stored and handled properly.

Warranty

Finnzymes Oy warrants that its products will meet the specifications stated on the technical data section of the data sheets, and Finnzymes Oy agrees to replace the products free of charge if the products do not conform to the specifications. Notice for replacement must be given within 60 days of receipt. In consideration of the above commitments by Finnzymes Oy, the buyer agrees to and accepts the following conditions:

- That this warranty is in lieu of all other warranties, express or implied;
- That **ALL WARRANTIES OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE ARE HEREBY EXCLUDED AND WAIVED**;
- That the buyer's sole remedy shall be to obtain replacement of the product free of charge from Finnzymes Oy; and
- That this remedy is in lieu of all other remedies or claims for damages, consequential or otherwise, which the buyer may have against Finnzymes Oy.

Exclusive terms of sale

Finnzymes Oy does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorised officer of Finnzymes Oy. Prices are subject to change without notice.

Recommended guidelines for safe use of the products

Finnzymes Oy recommends that the buyer and other persons using the products follow the Guidelines for Research Involving Recombinant DNA Molecules (NIH guidelines) Federal Register, July 5, 1994 (59 FR 34496) and any amendments thereto. Finnzymes Oy disclaims any and all responsibility for any injury or damage which may be caused by the failure of the buyer or any other person to follow said guidelines.

Research use only

Since these products are intended for research purposes by qualified persons, the Environmental Protection Agency does not require us to supply Premanufacturing Notice.

Notice to user

The information presented here is accurate and reliable to the best of our knowledge and belief, but is not guaranteed to be so. Nothing herein is to be construed as recommending any practice or any product in violation of any patent or in violation of any law or regulation. It is the user's responsibility to determine for himself or herself the suitability of any material and/or procedure for a specific purpose and to adopt such safety precautions as may be necessary.

Notice to Purchaser of Phusion® DNA Polymerases

Limited license (proofreading DNA polymerases). The purchase price of this product includes a limited, non-transferable license under U.S. and foreign patents (5,500,363 and 5,352,778) owned by New England Biolabs, Inc. to use this product. No other license under these patents is conveyed expressly or by implication to the purchaser by the purchase of this product.

Notice to Purchaser

Limited license. The purchase price of this product includes a limited, non-transferable license under U.S. and foreign patents owned by BIO-RAD Laboratories, Inc., to use this product. No other license under these patents is conveyed expressly or by implication to the purchaser by the purchase of this product.

This product is sold under license from Affibody AB, Sweden.

The quality system of Finnzymes Oy is certified according to standard SFS-EN ISO9001:2008.

Phusion® is a trademark or a registered trademark of Finnzymes Oy.

DyNAzyme™ is a trademark of Finnzymes Oy.

Affibody® is a registered trademark of Affibody AB, Sweden.

Version 1.7, November 2009