



## Animal tissues tested with Phire<sup>®</sup> Animal Tissue Direct PCR Kit

June 2010

Mouse (ear, tail, embryo, bone, thymus, tooth, fat, kidney, brain, heart, liver, spleen, lung)

Cultured mouse cells (fibroblasts)

*Drosophila* (wing, whole insect)

Zebrafish (fin)

Whitefish (fin and scale)

*Caenorhabditis elegans* (whole worm)

Mosquito (whole insect)

Aphis (whole insect)

Catfish (fin)

*Agrilus viridis* (leg, body)

Bear (muscle)

Lynx (muscle)

Reindeer (muscle)

Roe deer (muscle)

Moose (muscle)

Bear (muscle)

Rabbit (muscle)

Wolverine (hair, bone, claw)

Bovine (hair)

Dog (hair)

Siberian jay (feather)

Woodpecker (feather)

Check the latest updates at

[www.finnzymes.com/directpcr](http://www.finnzymes.com/directpcr)